

# prisma VENT40

Invasive and non-invasive ventilation.

With prisma VENT40 you'll have mobility and convenience at every level:

- pressure range up to 40 hPa
- use of leakage circuit system (15 mm / 22 mm)
- with autoS/T mode (autoEPAP, autoF) and target volume (100 – 2,000 ml)
- includes helpful functions for COPD treatment: AirTrap Control, maximum setting for expiratory ramp and Trigger Lockout
- · integrated port for oxygen feed
- three ventilation programs stored in device memory
- complete alarm management
- rechargeable battery with 10-hour capacity
- easy to use thanks to intuitive menus, simple operation and quick access

Quiet as a whisper – for the comfort and well-being of patients and partners.

### **Our Software Solutions**



LÖWENSTEIN

prisma CLOUD

The telemedical connection for experts and patients.
Information updated daily on therapy, remote setting of therapy parameters. Fully automated, from anywhere.

prismaTS

STATISTICAL PROPERTY OF THE PARTY OF THE PAR

Therapy software for experts. For detailed data analysis and remote settings.



## prisma VENT40

0.5 s - 4 s 0.2 s - 4 s auto (only Ti timed)

#### **Accessories**

<b>Humidifier</b> prisma VENT AQUA prismaAQUA	WM 100506 WM 29490
Leakage circuit 22 mm Ø, black 22 mm Ø, autoclavable 15 mm Ø	WM 23962 WM 24667 WM 29988

You'll find the complete line of accessories in the prisma VENT accessories catalog.

#### prismaHYBERNITE

breathing tube, heated

22 mm Ø WM 29067 15 mm Ø WM 29083

Inverter

12 V WM 24616

#### **Technical specifications**

Product class as per MDR (EU) 2017/745	IIa	Ti/Ti max Ti min, Ti max, Ti timed
Dimensions in cm (W x H x D):	21.8 x 17.5 x 21.8	Respiratory rate
Weight	2.4 kg	Target volume
Weight of internal rechargeable battery (for variant with battery)	0.63 kg	Trigger levels - Inspiration
Maximum air flow at 20 hPa	> 220 l/min	- Expiration
Operating time of internal rechargeable battery with the following settings: T mode, f = 20 /min, Ti = 1 s, PEEP = 4 hPA, Vt = 800 ml	> 10 hours	Pressure increase speed
Passive lungs: Resistance R = 5 hPa (l/s); Compliance C = 50 ml/hPa		Pressure decrease speed
IPAP pressure range	4 – 40 hPa ±1.2 hPa (±8% of setting)	
PEEP pressure range	4 – 25 hPa ±1.2 hPa (±8% of setting)	Maximum allowed flow with oxygen supply

	auto (orny 11 timea)
Respiratory rate	0 – 60 breaths/min
Target volume	100 ml – 2000 ml
Trigger levels	
- Inspiration	1 (high sensitivity) to 8 (low sensitivity)
- Expiration	95% to 5% of maximum flow in steps of 5%
Pressure increase speed	Level 1: 100 hPa/s
	Level 2: 80 hPa/s Level 3: 50 hPa/s
	Level 4: 20 hPa/s
Pressure decrease speed	Level 1: 100 hPa/s
	Level 2: 80 hPa/s
	Level 3: 50 hPa/s
	Level 4: 20 hPa/s
	max: slow pressure adjust- ment
Maximum allowed flow with oxygen supply	15 l/min

#### **Device designation**

#### Article no.

prisma VENT40 WM 29500-1110 prisma VENT40 with integrated rechargeable battery WM 29550-1110

Please refer to the currently valid price list.

CE 0197

**AUTHORISED IMPORTER** LAOS MEDITEK A-1, PARMESHWAR BUNGLOW, B/H ANNAPURNA RESTAURANT, JASHODANAGAR CR RD, AHMEDABAD-382445 PHONE: 9879512469





Löwenstein Medical Technology Kronsaalsweg 40 22525 Hamburg, Germany T: +49 40 54702-0 F: +49 40 54702-461 info@loewensteinmedical.de loewensteinmedical.com

#### Sales + Service

Löwenstein Medical Arzbacher Straße 80 56130 Bad Ems, Germany T. +49 2603 9600-0 F. +49 2603 9600-50 info@loewensteinmedical.com loewensteinmedical.com









